

# LEADTOOLS Medical Web Viewer

## Requirements

### Client Side:

- Google Chrome version 71.0 or later
- Microsoft Internet Explorer 11
- Microsoft Edge 14.95 or later
- Mozilla Firefox 65.0 or later
- Microsoft .NET Framework 4.0

### Server Side:

#### Minimum Requirements: Medical Viewer only (without 3D viewing)

- OS: Windows Server 2008 R2 or later  
Windows 7 SP1 or later  
Microsoft .NET Framework 4.5
- SQL: SQL Server 2008 R2 or later
- CPU: Intel Core 2 Duo E6600 or higher
- RAM: 4 GB (DDR3)
- Disk: 100 GB
- Graphics Card: N/A
- GPU Memory: N/A

#### Minimum Requirements: Medical Viewer with 3D (Software Rendering)

- OS: Windows Server 2008 R2 or later  
Windows 7 SP1 or later  
Microsoft .NET Framework 4.5
- SQL: SQL Server 2008 R2 or later
- CPU: Intel i3-2100 3.1GHz or higher
- RAM: 8 GB (DDR3)
- Disk: 500 GB
- Graphics Card: N/A
- GPU Memory: N/A

#### Minimum Requirements: Medical Viewer with 3D (Hardware Rendering)

- OS: Windows Server 2008 R2 or later  
Windows 7 SP1 or later  
Microsoft .NET Framework 4.5
- SQL: SQL Server 2008 R2 or later
- CPU: Intel i3-2100 3.1GHz or higher
- RAM: 8 GB (DDR3)
- Disk: 500 GB
- Graphics Card: DirectX 9 Compatible
- GPU Memory: 1 GB

### Recommended Requirements: Medical Viewer with 3D

- OS: Windows Server 2012 R2 or later  
Windows 7 SP1 or later  
Microsoft .NET Framework 4.5
- SQL: SQL Server 2008 R2 or later
- CPU: Intel XEON E5-2650v4 or higher
- RAM: 64 GB (DDR4)
- Disk: 1 TB SSD
- Graphics Card: DirectX 9 Compatible
- GPU Memory: 8 GB GDDR5 or higher